The Agilent 2100 Bioanalyzer

RNA Integrity Number (RIN) A standardized approach for RNA integrity assessment



qPCR Symposium Leipzig, 2005



Marc Valer Agilent Technologies

RNA LabChip kits Analysis of Total RNA Integrity





Current total RNA integrity measures

- UV offers no information on size distribution, and does "see" interferences of DNA, solvents and other contaminants. Sophisticated de-convolution is possible to individually quantity the components.
- Gels, have a limited linear range and ribosomal band intensity ratio had been used to differentiate good to bad samples. This measure is very subjective, dependent on the
- RNA degradation is a gradual process.
- Results have to be interpreted by visual inspection.
- Overlay of electropherograms only works well for samples with the same concentration.
- Instrument dependency in signal height



mRNA QC is addressed with smear



Problem Description



- The ratio of ribosomal bands is not sufficient to describe RNA integrity!
- RNA degradation is a gradual process.
- Results have to be interpreted by visual inspection.
- Overlay of electropherograms only works well for samples with the same concentration.
- Instrument dependency in signal height



Goal

- Find a tool that describes the integrity of a RNA sample better than the ribosomal ratio alone could do.
- Have a bioanalyzer result that describes the sample independently of individual interpretation.
- Characterize a RNA sample reproducibly independent on the instrument is was analyzed on and independently of who is performing the analysis
- Allow comparison of validation results across different laboratories



Introducing the RNA Integrity Number (RIN)

- The RIN is a software algorithm that allows classification of RNA integrity
- It extracts a number of characteristic features from the bioanalyzer electropherogram
- An adaptive learning process is used to "teach" the algorithm about the relative importance of the extracted features (based on many example electropherograms).
- The RIN returns a value that is characteristic for the integrity of a specific sample.



Approach – Flow chart





Data Preparation

- Initially limited to Eukaryote Total RNA Nano Assay, expanded to Eukaryote RNA pico Assay
- Input data: Random selection of mostly <u>human</u>, <u>mouse</u> and <u>rat</u> total RNA
- >1400 samples analyzed and classified by application specialists, origin of samples: microarray service laboratory and Agilent Labs
- Manual generation of quality labels, q = targets for learning
- Manual generation of anomaly labels, a (Ghost Peaks, wavy Baseline)
- Quality is a continuous variable, no separate classes by nature
 → Introduction of 10 discrete classes (1 to 10 in ascending quality)



Feature Generation

Approach: Training of regression models for peak positions and Intensities based on peak labels



Selection of the best feature combination



RIN Visualization







RNA Area:		7,007.2							
RNA Concentration:		2,248	ng/µl						
rRNA Ratio [28s / 18s]: 0.0									
RNA Integrity Number (RIN: 2.5									
<u>R</u> esults	<u>P</u> eak Table	Frag <u>m</u> ent Tal	ble <u>L</u> egend	<u>E</u> rrors					



RIN - Limitations

• What the RIN can do:

- Obtain an assessment of RNA integrity.
- Directly compare RNA samples
- Ensure repeatability of experiments

• What it CANNOT do:

- Predict the outcome of an experiment if no prior validation was done



RIN Visualization

	RNA Area: RNA Concentration:	26.3 75	ng/µ	I		Anomaly threshold settings
	rRNA Ratio [28s / 18s]:	1.9				🖃 : RNA Integrity Number
2	RNA Integrity Number (RIN):	96				Pre Region Anomaly Thres 0.5
	rant mogney framber (rany					55 Region Anomaly Thresh 0.3
						Fast Region Anomaly Thre 0.77
						Inter Region Anomaly Thre 0.21
						Precursor Region Anomaly 0.3
						Post Region Anomaly Thre 0.29
						Baseline Anomaly Threshold 0.56
						Ribosomal Ratio Anomaly T 0.89
						Unknown Sample Type Thr 0.83
						Marker Anomaly Threshold 0.88
	Results Peak Table Fi	rag <u>m</u> ent Tab	le	Legend	Errors	



Results



RIN Application – Assessment of RNA Integrity



RIN Application – Assessment of RNA Integrity



RIN Application – Directly Compare Samples same sample on different instruments



When testing an identical RNA sample on various instruments, identical RINs are obtained – within narrow limits 36 samples 3 instruments CV RIN: 3 % CV ribosomal ratio: 12 %



RIN Application – Directly Compare Samples same sample in different dilutions



When testing an identical RNA sample in various dilutions, identical RINs are obtained – within narrow limits 108 samples 3 dilutions CV RIN: 3 % CV ribosomal ratio: 17%

Agilent Technologies

RIN Application – Comparison to alternative integrity measures

When applying the degradometer tool to a random data set, a significant number of samples are "unclassified" (BLACK CODING)



The degradometer is based on the ribosomal ratio and the area corresponding to our "fast region". The degradometer classifies samples in 4 categories, when RIN data is clustered in blocks the data correlates.

Agilent Technologies

RIN Validation for Gene Expression Experiments





RNA QC in Routine Gene Expression Workflow



Continue with downstream experiment (microarray, real-time PCR, etc.)



The effect of RNA degradation





Labeled cRNA profile



RIN to DIN? aCGH sample QC amplification

Genomic DNA amplification with phi29 titration series (50ng, 0.05ng, 0.005ng & no input) analysez with the RNA 6000 Nano kit. The increasing non-specific peaks can be seen in the 30-40 second range and the decreasing specific peaks can be seen migrating in the 50-60 second range.





The same samples were analyzed on a 0.8% agarose gel and visualized with EtBr.

600ng of each sample was loaded

Lanes: 1-1kB marker, 2-50ng, 3-5ng, 4-0.5ng, 5-0.05ng, 6-0.005ng, 7-0.0005ng, 8-no input

Can not distinguish between the specific and non-specific amplification peaks.



Conclusion

- A software algorithm was presented that allows standardization of RNA integrity assessment
- Ribosomal ratio is important parameter for RNA integrity but many other factors have to be taken into account
- Initial results indicate that the software algorithm is independent of instrumentation and can ensure repeatability of results (e.g. RT-PCR or microarray)
- Availability of RIN-tool:
 - included in beta-SW version: downloadable via Agilent Web-Page
 - included in next official SW release: May 2005



Acknowledgment

- Software
 - Andreas Schroeder
 - Michael Leiber
 - Susanne Stocker
 - Thomas Ragg



- Application
 - Samar Lightfoot
 - Ruediger Salowsky
 - Marcus Gassmann
 - Christine Miller*
 - Rainer Wittig⁺

We would like extend special thanks to our collaboration partners Ambion Inc. and the German cancer research center as well as to Quantiom Bioinformatics. We would also like to thank all researchers who beta-tested the RIN and provided valuable

Agilent Technologies